

COVID-19 home test



COMPATIBLE DEVICE LIST

IVD No test component to be used inside the body except the Nasal Swab as directed.

www.ellumecovidtest.com

What phones are compatible with the Ellume COVID-19 Home Test:

The Ellume COVID-19 Home Test is compatible with most Bluetooth enabled smartphones.



Apple iPhones:

Minimum supported capabilities: an ARM64 iPhone device with Bluetooth running iOS 12 or later.

Inclusions:

- iPhone 5s
- iPhone 6
- iPhone 6 Plus
- iPhone 6s
- iPhone 6s Plus
- iPhone SE
- iPhone 7
- iPhone 7 Plus
- iPhone 8
- iPhone 8 Plus
- iPhone X
- iPhone XS
- iPhone XS Max
- iPhone XR
- iPhone 11
- iPhone 11 Pro
- iPhone 11 Pro Max
- iPhone SE (2nd Generation)
- iPhone 12
- iPhone 12 mini
- iPhone 12 Pro
- iPhone 12 Pro Max

The test is not compatible with iPads, iPod Touch, iPhone 5 (ArmV7s devices).

Android phones:

Minimum supported capabilities: an Android smartphone with Bluetooth 4.0 or later and XHDPI (720x1280px) resolution supporting Android API 23 (version 6.0) or later.

The test is not compatible with Android tablets.