

# COVID-19 home test



## COMPATIBLE DEVICE LIST

**IVD** No test component to be used inside the body except the Nasal Swab as directed.

[www.ellumecovidtest.com](http://www.ellumecovidtest.com)

**The Ellume COVID-19 Home Test is compatible with most Bluetooth enabled smartphones.**

### Apple iPhones:

Models released from 2013 onwards. The phone needs to be using an operation system iOS 10 onwards.

#### Compatible phones

(list up to date as of November 2020)

Models not on this list but released from 2013 and operation iOS 10 will be compatible.

- iPhone 5s
- iPhone 6
- iPhone 6 Plus
- iPhone 6s
- iPhone 6s Plus
- iPhone SE
- iPhone 7
- iPhone 7 Plus
- iPhone 8
- iPhone 8 Plus
- iPhone X
- iPhone XS
- iPhone XS Max
- iPhone XR
- iPhone 11
- iPhone 11 Pro
- iPhone 11 Pro Max
- iPhone SE (2nd Generation)
- iPhone 12
- iPhone 12 mini
- iPhone 12 Pro
- iPhone 12 Pro Max

**The Ellume COVID-19 Home Test is not compatible with iPads, iPod Touch, iPhone 5 (Armv7s devices).**

### Android phones:

Android smartphones (for example Samsung, Google, Motorola, Huawei, LG, Nokia, Alcatel, ZTE) with Bluetooth 4.0 or later. The phone needs to have XHDPI (720x1280px) resolution supporting Android API22 (version 5.1) or later.

**The Ellume COVID-19 Home Test is not compatible with Windows smartphones or Android tablets.**